

S-200/S-220/S-300 Spectrophotometer

INSTRUCTION MANUAL

Version V1.1



S-200/S-220/S-300 Spectrophotometer



Thank you for your purchase of S-200/S-220/S-300 spectrophotometer. In order to ensure proper use and your safety, please read this manual carefully and keep it well before using the instrument.

Information contained in this manual is subject to change without notice for product appearance and technical data. Enquiries are welcome.



IMPORTANT

Precautions on Electromagnetic Wave Interference

(1) Possible Electromagnetic Wave Interference Caused by This Instrument

When this instrument is used in a residential area or an adjacent area thereto, it may cause interference to radio and television reception.

To prevent this, use the specified system connection cables in strict accordance with the instruction manual.

The instrument is designed to minimize possible electromagnetic wave interference caused by it if the specified cables are connected properly.

However, there is no guarantee that electromagnetic wave interference will not be caused by the instrument.

If the instrument does cause interference to radio or television reception, which can be determined by turning off and on the instrument, the user is encouraged to try to correct the interference by one or more of the following measures:

- Increase separation between the instrument and radio/TV receiver.
- Connect the instrument to an outlet on a circuit different from that to which the radio/TV receiver is connected.

(2) Possible Electromagnetic Wave Interference Affecting This Instrument

If this instrument is used near an intense electromagnetic source, interference noise may be given to the instrument to incur an adverse effect on its performance or functionality.

To prevent this, use the specified system connection cables in strict accordance with the instruction manual.

The instrument is designed to minimize possible electromagnetic wave interference affecting it if the specified cables are connected properly.

However, there is no guarantee that electromagnetic wave interference will not occur in this instrument.

If the instrument does incur electromagnetic wave interference, which can be determined by turning on and off possible sources of electromagnetic wave interference nearby, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient the instrument.
- Increase separation between the instrument and possible sources of electromagnetic wave interference.

• Increase separation between the power cable of the instrument and possible sources of electromagnetic wave interference.

• Connect the instrument to an outlet on a circuit different from that to which possible sources of electromagnetic wave interference are connected.

• Confirm that any other device connected with the instrument is not



affected by electromagnetic wave interference.

Warranty on Product

The Model S-200/S-220/S-300 spectrophotometer is warranted to operate according to the specifications given in the instruction manual, provided it is used in accordance with the instructions described in the manual.

(1) Scope of Warranty

(a) Any parts which prove to be defective in design or workmanship during the warranty period will be repaired without charge.

(b) A substitute part may be used for repair, or replacement with an equivalent product may be made instead of repair.

(c) Such system components as a personal computer and printer to be updated frequently for improvement may not be available in original versions at the time of replacement.

(2) Warranty Period

One year from the date of initial installation.

(3) Availability of Technical Support Service

Technical support service for this instrument is available within regular working hours on workdays predetermined by us.

(4) Limitations and Exclusions on Warranty

Note that this warranty is void in the following cases, even if they occur within the warranty period.

(a) Failure due to operation at a place not meeting the installation requirements specified by us

(b) Failure due to power supply voltage/frequency other than specified by us or due to abnormality in power supply

(c) Corrosion or deterioration of the tubing due to impurities contained in reagent, gas, air or cooling water supplied by the user

(d) Corrosion of the electric circuits or deterioration of the optical elements due to highly corrosive atmospheric gas

(e) Failure due to use of hardware, software or spare parts other than supplied by us

(f) Failure due to improper handling or maintenance by the user

(g) Failure due to maintenance or repair by a service agent not approved or authorized by us

(h) After disposal of this instrument, or after its resale without our approval

- (i) Failure due to relocation or transport after initial installation
- (j) Failure due to disassembly, modification or relocation not approved by us

(k) Consumables, and failure of parts that have reached the end of specified useful life

(I) Failure of parts excluded from the warranty in the instruction manual or other documents

(m) Failure due to acts of God, including fire, earthquake, storm, lightning, social disturbance, riot, crime, insurrection, war (declared or undeclared), radioactive pollution, contamination with harmful substances, etc.

(n) Failure of the hardware, or damage to the system software, application software, data or hard disk due to computer virus infection

(o) Failure of the personal computer connected with the instrument, or damage to the system software, application software, data or hard disk due to power interruption or momentary power voltage drop caused by lightning or the like

(p) Failure of the personal computer connected with the instrument, or damage to the system software, application software, data or hard disk due to disconnection of main power to the personal computer without taking the specified normal shutdown procedure

(5) Disclaimer of Warranty

(a) Any express warranties other than the explicit conditions indicated in (1) are excluded from this warranty.

(b) Any other implied warranties of merchantability and fitness for a particular purpose are not included in this warranty. No liability is assumed for direct or indirect damages arising out of explicit or implied warranties.

(c) Oral or written information or advice given by our dealers, distributors, agents or employees without our express permission shall not create a warranty or in any way increase the scope of this warranty.

Installation, Relocation and After-sale Technical Service

Installation of this instrument shall be carried out by or under supervision of qualified service personnel of Boeckel + Co (GmbH + Co) KG or its authorized service agent.

Before installation of the instrument, the user shall make preparations for satisfying the installation requirements in accordance with the instruction manual.

If relocation of the instrument becomes necessary after initial installation (delivery), please notify your local sales representative or nearest service office of Boeckel + Co (GmbH + Co) KG.

Disposal this instrument

When you discard equipment, please check and discard a related statute etc. or ask the service section of Boeckel + Co (GmbH + Co) KG.

Other Precautions

(1) Handling of Chemicals and Samples

(a)The user is responsible for following relevant legal standards and regulations in handling, storage and discarding of chemicals and samples used in analytical operations of this instrument.

(b) Reagents, standard solutions and accuracy-control samples shall be handled, stored and discarded as instructed by the respective suppliers.



(2) Notice on Instruction Manuals

(a)Information contained in the instruction manuals furnished with the instrument is subjected to change without notice for product improvement.(b) This manual is copyrighted by Boeckel + Co (GmbH + Co) KG with all rights reserved.

(c) No part of this manual may be reproduced or transmitted in any form or by any means without our express written permission.





CAUTION

For your safety please read the following precautions carefully before using the instrument.

General Safety Guidelines

- For safe handling of this product, please follow the instruction procedure in the manual for this product.
- Pay special attention to follow all the hazard warnings on the product and in the manual. Failure to do so can cause injury to you or damage to the product.
- After installation, please do not move the equipment. A vibration might affect the adjustment of the product.
- The hazard warnings, which appear on the warning labels on the product or in the manual, have one of the following alert headings consisting of an alert symbol and a signal word, DANGER, WARNING, or CAUTION.

DANGER:	Indicates an imminently hazardous situation that, if not				
	avoided, will result in death or serious injury. (It does not apply to this equipment.)				
	Indicates a potentially hazardous situation that, if not				
	avoided, can result in death or series injury.				
	Indicates a hazardous situation that, if not avoided, will or				
	can result in minor to moderate injury, or serious damage to the product.				
A: The alert s	symbol shown precedes every signal word for hazard				
warnings,	warnings, and appears in safety related descriptions in the manual.				

In addition, the following "Attention" and "Note" are not directly related to the safety of a person:



- \triangle Attention: It is used to present warnings, which are not directly related to personal injury hazards. It is used to indicate prevention against damage to the equipment.
- ◇ Note: This is used to indicate instructions that enable you to operate the equipment accurately and perform accurate measurements.

General Safety Guidelines (Continued)

Before using

- Before using this product, please make sure you read and understand the instructions.
- Please keep this manual in a safe and easily accessible place so that you can use it when necessary.
- Please make sure to use this product properly and follow the instructions as specified in this manual.
- Please make sure to understand and follow the instructions regarding safety in this manual.
- If you do not follow the instructions in this manual, an inaccurate analysis may result or bodily injury may occur.
- Because of danger, please make sure not to modify or alter the product, make sure not to use unspecified parts, and make sure not to operate the equipment by removing/defeating the safety device(s).
- When using chemicals, please make sure to ventilate the room well. If there is not enough ventilation, it may be hazardous to your health.
- Although we have carefully considered the instructions written on the products and manuals, it is possible for an unexpected event to occur. When operating the equipment, aside from following the instructions, be very cautious.

General Safety Guidelines (Continued)

Precautions for Installation • Maintenance • Relocation and After Sale Technical Service

- Before installation, confirm that there are no missing items or standard accessories. If there is something missing or damaged, or you have noticed any problems, please contact our nearest representative.
- Operating the equipment without a standard part can damage the equipment and cause safety concerns. If that occurs, please follow the instruction of the installer.
- Installation of this instrument shall be carried out by or under supervisions o personnel of Boeckel + Co (GmbH + Co) KG or its authorized service agent
- When relocation of this instrument becomes necessary after initial installation (delivery), please notify your local Boeckel + Co (GmbH + Co) KG sales representative or nearest Boeckel + Co (GmbH + Co) KG service office. Technical support service for this instrument is available from service agent approved or authorized by Boeckel + Co (GmbH + Co) KG within regular working hours or workdays.
- Please do not perform any other operations that are not included in the manual. If any problem occurs with the equipment, please contact the agent from whom you purchased it or the service department of Boeckel + Co (GmbH + Co) KG.



WARNING: Poisoning from Organic Solvent Gas

Handling Organic Solvents

■ The organic solvent vapor may be harmful to your health.



WARNING: Eye Injury from Organic Solvents

Handling Organic Solvents

Please wear protective glasses when using organic solvents. If the organic solvent should get into your eye, flush your eye immediately under running water for at least 5 minutes while keeping your eyelids open. See a physician for appropriate treatment.

WARNING: Electrical Shocks from Improper Grounding

When wiring the personal computer, power supply for the thermostatic cell holder and the like, please make sure to use the 3-prong wire (with ground) provided.



WARNING: Electrical Shock from Contact

High voltage is used inside of the equipment, be sure to turn off the unit before connecting the power cord.

CAUTION: Burns from high Temperature

- The lamp will become very hot during operation.
- Make sure that the instrument is switch off, the power cord is pulled off and the D2 lamp and the tungsten lamp is cooled off when replacing the lamp.

CAUTION: Fatigue due to Prolonged Work

- Viewing the display in your work can cause eye and physical fatigue if you continue to work in the same posture for extended periods.
- When working with the display for a prolonged period, for your health, make sure to take breaks for 10 to 15 minutes every hour in order to rest your eyes and body.

CAUTION: Indoor ventilation

- If using UV zero detection system in a small room for a long time, it may cause indoor nitrogen concentration increases and the oxygen concentration drops, which will affect human health.
- Please install exhaust fan or open the windows frequently to maintain good indoor ventilation.



Electricity

- (a) The voltage for the Spectrophotometer system and personal computer must be a single-phase AC 100V to 240V; Variations in the voltage and noise generated in the power line will cause adverse effects on the spectrophotometer and may also cause accidents.
- (b) Please make sure that grounding is provided together with the power supply wires, and make sure that it is connected with a grounding resistance of less than 4Ω . Defective grounding may not only cause lower resistance against noise from the outside but it can also cause the Mass Spectrometer and personal computer to generate static electricity, which may involve the danger of electrical shocks.
- (c) A high voltage circuit is used inside the Spectrophotometer. Do not open the covers when this circuit is operational because of the danger from electrical shock.

Fire Extinguishers

(a) Do not smoke or use fire within 3 meters of the Spectrophotometer

(b) Make sure to keep a fire extinguisher near the Spectrophotometer at all times. Obtain an ABC Powder extinguisher that can be used for normal fires, oil fires, and electrical fires.



Functional characteristics

This series includes 2 models. S-200 is visible spectrophotometer, S-220/S-300 is UV-VIS spectrophotometer.

Feature	
Small and beautiful appearance	4.3 inches, 480×272 resolution, TFT colorful liquid crystal touch screen size: like A3 paper, 400(W)×280(D)×160(H) mm weight: 4 kg
Easy operation	One-button operation interface into the measurement function High-speed wavelength move, arrive to any specified wavelength within 1 second For S-200 spectrophotometer, when sampling interval is 0.2nm, its scanning speed is 2400nm/min. Unique application method manager changes the original function-oriented operation method and user sample measurement method operation mode.
Energy saving and environmental protection	High conversion efficiency switching power supply, 100-240V AC voltage input, and complete riddance of low grid voltage S-200 spectrophotometer uses 2000 hours lifetime imported halogen lamp S-220/S-300 spectrophotometer uses low-power, high brightness and 10 ⁹ times lifetime pulses xenon lamp
More ports	Serial printer port for thermal printer USB port for PC SD card port saves data and measurement methods Options port for connecting and controlling several options
More options	Auto 5-cell holder Auto sample sipper Flow cuvette holder 10, 20, 30, 50, 100mm rectangular long-path cuvette holder Micro-cuvette holder Tube holder (only for S-200) Electronic thermostat TC cuvette holder (only for S-220/S-300)



Instructions for installation and use

The instrument



Keypad



The instrument uses touch display screen, all operation can be by touch screen. Two shortcut key further enhance easy operation.



BACK. Return to the previous menu.

MEASURE. Begin to a new measurement.



Operation environment

Power	Voltage: AC 100V-240V, Frequency: 50/60Hz±1Hz, Capacity: more than 200W Grounding line resistance of 10Ω or less is required The input power mutations, without interference from other large electrical equipment				
Operating temperature	5-35℃				
	In order to perform a measurement under the most stable condition, we recommend				
	that the instrument is used in an air conditioned room of 20-25°C.				
	Storage temperature:-20-55°C				
Operating	Less than 85%.				
humidity	We recommend that the instrument is used under $30\% \sim 70\%$ humidity.				
	Storage humidity: less than 85%。				
Workbench	Width: more than 550mm, Depth: more than 400mm, Load bearing: more than 15kg, horizontal and reliable workbench.				
	If there are other auxiliary devices, it is necessary to enlarge workbench area. More than 200mm space at the both sides of the instrument, avoid being close to the wall.				
	Height of the instrument: about 160mm. In order to easy operation, please choose suitable workbench.				
Connect with	Connect with PC (additional optional software and computer). Indirect specify type				
other equipment	thermal printer (options) and other accessories.				
	When main unit connects with above equipments, please switch off the power of the main unit and other equipments.				
	All grounding cords are in good condition and can be connected with the grounding line of the main unit.				



Specification (I)

S-200					
Wavelength range	320-1100nm				
Wavelength controlled	0.2nm				
variable					
Wavelength accuracy	±1nm				
Wavelength repeatability	≤0.5nm				
Transmittance accuracy	±0.5%T (NIST 930 Filter)				
Transmittance repeatability	0.2%T				
Baseline flatness	±0.002Abs (330-1090nm)				
Noise level	≤0.001Abs (500nm)				
Baseline stability	≤0.001Abs/h (500nm, 2 hours after warm up)				
Spectral bandwidth	6nm±1.2nm				
Stray light	≤0.5%T				
Wavelength scan speed	2400nm/min (0.2nm sampling interval, without filter				
	switchover)				
Wavelength move speed	To any specified position within 1 second				
Light source	WI lamp				
Detector	silicon photodiode				
Display screen	4.3 inches, 480×272 colorful touch LCD screen				
Printer	specified 80-column thermal printer (series port)				
Metering mode	Single beam				
Dimension (W×D×H)	400×280×160mm				
Weight	About 4 kg				
Power requirement	AC, 100-240V, (50/60Hz)				
Power consumption	100VA				
Communication ports	Serial printer port connects thermal printer				
	USB port connects PC				
	SD card port saves data and measurement methods				
	Accessories port connects and controls several options				
Optional Accessories	Auto 5-cell holder				
	Auto sample sipper				
	Flow cuvette holder				
	10, 20, 30, 50, 100mm long-path cuvette holder				
	Micro-cuvette holder				
	Tube holder				



Specification (II)

S-220/S-300					
Wavelength range	190-1000nm				
Wavelength controlled	0.2nm				
variable					
Wavelength accuracy	±2nm				
Wavelength repeatability	≤1nm				
Transmittance accuracy	±1%T (NIST 930 Filter)				
Transmittance repeatability	0.5%T				
Baseline flatness	±0.005Abs (200-990nm)				
Noise level	≤0.005Abs (250nm)				
Baseline stability	≤0.005Abs/h (250nm, after 2 hours warm up)				
Spectral bandwidth	5nm±1nm				
Stray light	≤0.5%T				
Wavelength scan speed	300nm/min (0.2nm sampling interval, without filter				
	switchover)				
Wavelength move speed	To any specified position within 1 second				
Light source	pulsed-xenon lamp				
Detector	silicon photodiode				
Display screen	4.3 inches, 480×272 colorful touch LCD screen				
Printer	specified 80-column thermal printer (series port)				
Metering mode	Single beam				
Dimension (W×D×H)	400×280×160mm				
Weight	About 4 kg				
Power requirement	AC, 100-240V, (50/60Hz)				
Power consumption	100VA				
Communication ports	Serial printer port connects thermal printer				
	USB port connects PC				
	SD card port saves data and measurement methods				
	Accessories port connects and controls several options				
Optional Accessories	Auto 5-cell holder				
	Auto sample sipper				
	Flow cuvette holder				
	10, 20, 30, 50, 100mm long-path cuvette holder				
	Micro-cuvette holder				
	Electronic thermostat TC cuvette holder				



The operation interface



Shortcut appear on every interface and can quickly switch to different measurement functions.

Measurement appear on every measure interface for basic measurement and setting

Data Process appear on scan interface. For the results of the scan data processing and according to different scanning characteristic, automatically vary corresponding processing ways.

Graph display buttons appear only in scan interface for graphical display changes

Status bar show in every interface to display current time, wavelength and readings in transmission or absorbance

Clicking the wavelength in the Status bar, it will go to the setting screen for wavelength setup.

Clicking the **reading** in the **Status bar**, it will go to the large display measurement interface to achieve simple transmittance and absorbance read function.

The value in the status bar of S-200 is real-time reading value of current wavelength.

The value in the status bar of S-220/S-300 is the last time measured value.



Word input keyboard

HOME			500.0nm	0.000ABS	12/06/28 09:35
ABS/%T	DATA	MODE P	ATHLENGTH	SAMPLE NAME	UL NUMBER
WL SCAN	ABS		10	5-01	3
TIME SCAN	0 1	T e T	рт	V II	I O P
CONC	A	S D	FG	HJ	KL
MANAGER	•	z x	C V	BN	MC
SETTING	123		Ĺ		+

This keyboard appears when there is a need for word input, for example, sample name, user and document name.



Numeric input

HOME		500.0n	m 0.000	ABS 12/0	6/28 09:35
ABS/%T					
WL SCAN		•	С		
TIME SCAN	0.0	1	2	3	
CONC		4	5	6	
MANAGER		7	8	9	
SETTING	₹T	•	0	_	

Numeric input appears at the area where can input digit, for example, wavelength value and scanning time.



Main Menu



F	unction	of	each	key

ABS/%T:	Photometric (single/ multiple wavelength)
WL SCAN:	Wavelength scan
TIME SCAN:	Time scan
CONC:	Concentration measurement
MANAGER:	File manger - Saved data & Saved method
SETTING:	System Setting



The large display can be used to measure single wavelength's absorbance or transmission. Click the value in the status bar in any screen can enter the large character measure mode.

ZERO: Blank sample

READ: Start a new measurement (only for S-220/S-300 spectrophotometer)

PRINT: Print results



Click wavelength in the status bar, input wavelength, the instrument will go to the select wavelength for measurement. (GOTO λ function)



Single/Multiple wavelength measurement



Photometry function can measure up to o selected wavelength s absorbance or transmission. Each measurement file can store up to 200 results.

READ:Start measurementSET:Parameter settingZERO:Blank sampleLOAD:Load data files or method filesPRN:Print results



Single/Multiple wavelength measurement setting

HOME	 Ø 		500.0nm	0.000A	BS 1	.2/06/28	09:35
ABS/%T	DATA MODE	PATH 10	ILENGTH	SAMPLE S-01	NAME	WL NUM	BER
WL SCAN							
	WL1(nm)		WL2 (:	nm)		WL3(nm)	
TIME SCAN	440	440		546		635	
CONC							
MANAGER				_			
SETTING		OK			Save		

DATA MODE:	%T or ABS mode
PATHLENGTH:	Cuvette pathlength
SAMPLE NAME:	Sample name input (less than 8 characters)
WL NUMBER:	Wavelength number (up to 6)
WL1~WL6:	Measured wavelength selection
OK:	Confirm setting
Save:	Save parameter



Wavelength Scan



Wavelength scan allow scanning the sample transmission, absorbance and energy (for service purpose).

READ:	Start measurement
SET:	Parameter setting
ZERO:	Blank sample
LOAD:	Load data files or method files
PRN:	Print results



Wavelength scan parameter setting

HOME		500.0nm 0.000AB	\$ 12/06/28 09:35
ABS/%T	DATA MODE	SAMPLE NAME	PEAK THRESHOL
	ABS	S-01	0.01
WL SCAN			
	START WL (nm)	END WL (nm)	INTERVAL (nm)
TIME SCAN	500	1000	AUTO
CONC	Y-MAX	Y-MIN	PATHLENGTH
	0.01	-0.01	10
MANAGER			
SETTING	OK		lave

DATA MODE:	%T or ABS mode
SAMPLE NAME:	Sample name input (less than 8 characters)
PEAK THRESHOL:	Select Y-axis Peak threshold setting
START WL(nm):	Starting wavelength
END WL(nm):	End wavelength
INTERVAL(nm):	Scan interval. Select from AUTO, 0.2nm, 0.5nm, 1nm, 2nm, 5nm. AUTO- scan
	range < 500nm, time interval is 0.5nm. scan range > 500nm, time interval is 1nm.
Y-MAX:	Y-axis maximum value of scan graph
Y-MIN:	Y-axis minimum value of scan graph.
PATHLENGTH:	Cuvette pathlength
OK:	Confirm setting
Save:	Save parameter



Wavelength scan results and data processing



Data processing icons

Peak and valley value. The determination value is based on the set threshold.			
Smoothing function < 6 smoothness			
Calculation function including Addition, subtraction, multiplication and division			
Peak area calculation function. According to the setting wavelength range, calculate peak area			
Wavelength tracking. Display wavelength and value			
Graph auto zoom button.			
Zoom in			
Zoom out			



Time scan



Time scan can scan the sample transmission, absorbance with a specific wavelength and within set time.

READ:	Start measurement
SET:	Parameter setting
ZERO:	Blank sample
LOAD:	Load data files or method files
PRN:	Print results



Time scan parameter setting

HOME		500.0nm 0.000AE	\$ 12/06/28 09 : 35
ABS/%T	DATA MODE	s	AMPLE NAME
WL SCAN	ADV		~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~
	WL (nm)	SCAN TIME(s)	INTERVAL(s)
TIME SCAN	500	300	AUTO
CONC	Y-MAX	Y-MIN	PATHLENGTH
	0.01	-0.01	10
MANAGER			
	OK		Save
SETTING			

DATA MODE:	%T or ABS mode
SAMPLE NAME:	Sample name input (less than 8 characters)
WL(nm):	Scanning wavelength
SCAN TIME(s):	Scan time. (unit: s). more than 60s, max: 60000s
INTERVAL(s):	Scan interval. Select from AUTO, 0.5s, 1s, 2s, 5s, 10s, 100s. AUTO interval, when
	scan time is < 1000s, scan interval is 1s. When scan time is 1000-2000s, scan
	interval is 2s. When scan time is 5000-20000s, scan interval is 10s. When scan time
	is > 20000s, scan interval is 100s.
Y-MAX:	Y-axis maximum value of scan graph
Y-MIN:	Y-axis minimum value of scan graph.
PATHLENGTH:	Cuvette pathlength
OK:	Confirm setting
Save:	Save parameter



Time scan results and data processing



Data processing icons

Kinetics:	Dynamics calculation function. According to setting time range, make dynamics calculation
Smooth:	Smoothing function < 6 smoothness
Calcu:	Calculation function, including Addition, subtraction, multiplication and division
Area:	Peak area calculation function. According to the setting wavelength range, calculate peak area
h	Time tracking. Display time and value
e.	Graph auto zoom button.
Œ	Zoom in
e	Zoom out



Concentration measurement



Concentration or specified wavelength can be measured. Op to 5 wavelengths can be chosen. Up to 9 standard and 5 replicates can be recorded. Standard curve and its K-factor value are displayed. The standard curve data can be input for known concentration curve. Standard curve can be saved in the SD card. Each file can store up to 200 results.

READ:	Start measurement
SET:	Parameter setting
ZERO:	Blank sample
LOAD:	Load data files or method files
PRN:	Print results



Concentration measurement setting (I)

HOME		500.0nm	0.000ABS 1	.2/06/28 09:35
ADC /AT	DATHI FNCTH	SAMDIE	N & MF	
AD3/31	10	S-01	MALL	1
WL SCAN				
	WL1 (nm)			
TIME SCAN	440			
COMC		TEDATTON	STD MIMBER	DEDLICATES
CONC	ng/nl S	TD	3	1
MANAGER				
SETTING	ОК	Nex	t S	ave
SETTING				

PATHLENGTH:	Cuvette pathlength
SAMPLE NAME:	Sample name input
WLNUMBER	wavelength number (up to 3)
(nm):	
WL1~ WL3 (nm):	Measured wavelength selection
UNIT:	Sample's unit display
CALIBRATION:	Concentration calculation method. Select measurement with standard sample for
	fresh standard curve or input with standard curve K-factor.
STD NUMBER:	Number of Standard sample. Up to 9 sample.
REPLICATES:	Standard samples replicates number. Up to 5 replicates
OK:	Confirm setting
Next:	Next page
Save:	Save parameter



Concentration measurement setting (II)

HOME	I	500.0nm 0.00	00ABS 12/06/28	09:35
3DC /2T	(ITT) 1	em o	erm 2	
AD 3/31		10	20	
WL SCAN				
TIME SCAN				
CONC				
MANAGER	0			
			Back	
SETTING			D'den	

When selected measurement with standard sample for fresh standard curve- CALIBRATION , select STD, a setup page will be shown.

STD1~STD9: input the known calibration standard value

Concentration measurement setting (III)

HOME	 Image: A set of the set of the	500.0nm	0.000ABS	12/06/28	09:35
ABS/%T			v	.1	_
	0		1		
WL SCAN					
TIME SCAN					
CONC					
MANAGER				_	
SETTING		OK	Back		

When selected measurement with input standard curve K-factor- CALIBRATION , select K-factor, a setup page will be shown.

- K1: Slope value of standard curve
- K0: Offset value of standard curve
- OK: Confirm setting
- BACK: Previous page



File manager – Method (I)



Method Manager is designed for managing all method files. SD card must be inserted for storing and opening the files.

ABS/%T:	Saved Photometric methods folder
WL SCAN:	Saved Wavelength scan methods folder
TIME SCAN:	Saved Time scan methods folder
CONC:	Saved Concentration measurement methods folder
Open:	Open measure method
New:	Renew selected measure method



File Manager – Method (II)

HOME		500.0nm	0.000ABS 12/	06/28 09:35
DATA	2		Name: Method:	01 ABS/%T
METHOD	01 02 .PRP .PRF SSS .PRP	SA P.PRP	Data Mode: WL1: WL2: WL3:	ABS 440 546 635
	Open	Delete	Next	Back

After selected memory noncer in the previous menu. The previous saved individual memory in the corresponding method's folder

Open:	Open selected method file
Delete:	Delete selected method file
Next:	Next page
Back:	Previous page

Method file type

- *. PRP Photometric method file
- *. PRW Wavelength scan method file
- *. PRT Time scan method file
- *. PRC Concentration measurement method file



File Manager - Data



Data manager function is designed for managing all data files. SD card must be inserted for storing and opening the files.

- Open:Open selected method fileDelete:Delete selected method fileNext:Next page
- Back: Previous page

Data document type

- *. PHT Photometric data file
- *. WSC Wavelength scan data file
- *. TSC Time scan data file
- *. CON Concentration data file


System setting

HOME	V	500.0nm	0.000ABS	12/06/28	09:35
SYSTEM	USER	NAME	CON	IPANY	
TIME					
	HOME	SCREEN	INIT	IALIZE	
GLP/GMP	RTD				
	FACTO	ORY RESET			

USER NAME:	Input user name. (Less than 20 characters)		
COMPANY:	Input company name. (Less than 20 characters)		
HOME SCREEN:	Select the function showed in home page.		
	RTD: Zoom in photometric.		
	ABS/%T: Photometric function.		
	WL SCAN: Wavelength scan.		
	TIME SCAN: Time scan interface.		
INITIALZE:	System startup i.e. self checking. During the process, the sample chamber must keep closing.		
FACTORY RESET:	Restored to factory default settings. Sample name, operation name and company name in any interface will be empty.		



Time setting

HOME	 Image: A start of the start of	500.0nm	0.000ABS	12/06/28	09:35
SYSTEM	YEAR		М	ONTH	
TIME	12		0	6	
	DAY		Н	OUR	
GLP/GMP	28		9		
	MINUTE		F	ORMAT	
	35		1	TY/MM/DD	

There is a power-off storage function (~ 96 hours), the user may need to reset time after this period. Setting data did not affect.



GLP/GMP performance self-check

HOME		500.0nm	0.000ABS	12/06/28	09:35
SYSTEM	NOISE		STA	BILITY	
TIME					
GLP/GMP	BASELINE	FLATNESS	HAF	DWARE	
		Prii	at		

GLP/GMP is used to validate instrument's performance.

NOISE:	Equipment background noise level test
STABILITY:	Stability test, for best result need to preheat for 2 hours
BASELINE	Baseline flatness test
FLATNESS:	
HARDWARE:	Hardware self-check, including PCB, lamp
Print:	Print test result



S-300 interface

S-300 is only for analysis of nucleic acids and proteins. Besides the special function measure interface, it also own measure functions of traditional spectrophotometer.

S-300 mian interface



Click the icon, enter into corresponding measure interface.

DNA/RNAnucleic acids analysisPROTEINprotein determinationOD600Bacterial cell culture measurementSPECTRUMTraditional spectrophotometer functionMANAGERmeasure method manager including data-savedSETTINGparameter setting



S-300 nucleic acids analysis interface

HOME	500.0	nm 0.000ABS 12/06/28	09:35
DSDNA		UL (nm) OD	READ
SSDNA	CONCENTRATION	230 260	SET
RNA	CONCENTRATION	280 320	ZERO
OLIGO	ug/ml	A260/A280	LORD
			PNR

There are 4 methods in the DNA/RNA nucleic acids analysis interface, dsDNA, ssDNA, RNA, OLIGO. Click button to select measure method.

Measure buttons

- READ being to measure
- SET parameter setting
- ZERO zero absorbance (blank sample) at all setting wavelength point
- LOAD load original measure results or measure method document
- PRN print data



S-300 nucleic acids parameter setting interface

HOME	 Image: A start of the start of	500.0nm	0.000ABS	12/06/28	09:35
	EA CTOD		Pac	VCDOIND	
DSDNA	50		0	N	
SSDNA					
	PATHLENGTH		DI-	FACTOR	
RNA	10		1		
_			S 9 1	DIE MAME	
OLIGO	ua/m]		SAI S-	-01	
	OK		Sat	70	

FACTOR	regulate parameters
BACKGROUND	background calibration
PATHLENGTH	optical length of cuvette (optical length is less than 10mm)
DI-FACTOR	dilution factor
UNIT	concentration units
SAMPLE NAME	sample name
OK	Finish parameter setting, back to measure interface automatically
Save	save parameter (SD card)
PATHLENGTH DI-FACTOR UNIT SAMPLE NAME OK Save	optical length of cuvette (optical length is less than 10mm) dilution factor concentration units sample name Finish parameter setting, back to measure interface automaticall save parameter (SD card)



S-300 protein analysis interface



There are 5 analysis methods in protein determination interface, PROTEINUV, BRADFORD, BSA, LOWRY, BIURET. Click button to select measure methods.

PROTEINUV functions

Measure buttons

- READ begin to measure
- SET parameter setting
- ZERO zero absorbance (blank sample) at all setting wavelength point
- LOAD load original measure results or measure method document
- PRN print data

S-300 protein analysis interface PROTEINUV parameter setting interface

HOME	V	500.0nm 0.0004	ABS 12/06/28 09:35
PROTEINUU	A280 FAUTOR	A	260 FAUTUR
	1.45		0.76
BRADFORD			
	PATHLENGTH	BACKGROUND	PROTEIN
BSA	10	ON	BSA
LOWRY	DI-FACTOR	UNIT	SAMPLE NAME
	1	ng/nl	S-01
BIURET			
	07	_	Satto
	ON		

A280 FACTOR	correction factor at 280nm position
A260 FACTOR	correction factor at 260nm position
PATHLENGTH	optical length of cuvette (optical length is less than 10mm)
BACKGROUND	background calibration
PROTEIN	select measure methods
DI-FACTOR	dilution factor
UNIT	concentration units
SAMPLE NAME	sample name
OK	Finish parameter setting, back to measure interface automatically
Save	save parameter (SD card)



S-300 protein analysis interface

Select parameter of MOLAR EXT C (self-defining) protein measure method in PROTEINUV parameter setting interface

HOME	v	500.0nm 0.000ABS	12/06/28 09:35
PROTEINUU	UL (nm)	MOLAR EXT C	MOLECULAR W
	200	47790	69323.4
BRAD FORD			
BSA			
LOWRY			
BIURET			
	OK	Ba	ck

WL(nm)	Wavelength value
MOLAR EXT C	Molar extinction coefficient correction factor
MOLECULAR W	Molecular weight. According to the two correction factors, calculate correction factor at
	A280 FACTOR automatically
OK	Finish parameter setting, back to measure interface automatically
Back	Back to the prior parameter setting interface



S-300 protein analysis interface

Select parameter of EXT COEFF (self-defining) protein measure method in PROTEINUV parameter setting interface

HOME	V	500.0nm	0.000ABS	12/06/28	09:35
PROTEINUU	WL (nm)		EXT	COEFF	
	200		0.6	59	
BRADFURD					
BSA					
LOWRY					
BIURET					
	OK		Bac	k	

WL(nm)	Wavelength value
--------	------------------

- EXT COEFF Extinction coefficient correction factor. According to the correction factors, calculate correction factor at A280 FACTOR automatically
- OK Finish parameter setting, back to measure interface automatically
- Back Back to the prior parameter setting interface

S-300 protein analysis interface BRADFORD measure interface



Measure buttons

- READ begin to measure
- SET parameter setting
- ZERO zero absorbance (blank sample) at all setting wavelength point
- LOAD load original measure results or measure method document
- PRN print data

S-300 protein analysis interface BRADFORD parameter setting interface

HOME	I	500.0nm 0.000AB	s 12/06/28 09:35
PROTEINUV	UL (nm)	P.	ATHLENGTH
BRADFORD	595		10
Dimbrond	CALIBRATION	STD NUMBER	REPLICATES
BSA	STD	3	1
LOWRY	DI-FACTOR	UNIT	SAMPLE NAME
	1	ng/nl	S-01
BIURET	077		2000
			1008

WL(nm)	Wavelength value
PATHLENGTH	optical length of cuvette (optical length is less than 10mm)
CALIBRATION	concentration regression way
STD NUMBER	concentration standard sample number
REPLICATES	measure times
DI-FACTOR:	dilution factor
UNIT	concentration units
SAMPLE NAME	sample name
OK	Finish parameter setting, back to measure interface automatically
Save	save parameter (SD card)



S-300 protein analysis interface

CALIBRATION (concentration) standard sample setting interface in BRADFORD

parameter setting interface

HOME	V	500.0nm 0.000AB	\$ 12/06/28 09:35
PROTEINUV	STD1	STD2	STD3
BRADFORD	0		20
BSA			
LOURY			
BIURET			
	OK	В	ack

- STD1~3 standard sample concentration setting (According to the number of concentration in the upper layer interface, it will display corresponding setting)
- OK Finish parameter setting, back to measure interface automatically
- Save save parameter (SD card)



S-300 bacterial cell culture measurement interface OD600 measure interface



Measure buttons

- READ begin to measure
- SET parameter setting
- ZERO zero absorbance (blank sample) at all setting wavelength point
- LOAD load original measure results or measure method document
- PRN print data

S-300 bacterial cell culture measurement interface

OD600 parameter setting interface (I)

HOME	V	500.0nm	0.000ABS	12/06/28	09:35
0D600	WL (nm)		PAT	HLENGTH	
	600		10		
	UNIT		SAM	PLE NAME	
	OD		S-()1	
	OK		Sat	ле	
			14.5		

WL(nm)	Test wavelength
PATHLENGTH	optical length of cuvette (optical length is less than 10mm)
UNIT	concentration units
SAMPLE NAME	sample name
OK	Finish parameter setting, back to measure interface automatically
Save	save parameter (SD card)

S-300 bacterial cell culture measurement interface

OD600 parameter setting interface (II)

HOME	 Image: A start of the start of	500.0nm	0.000ABS	12/06/28	09:35
	1 A.29				
0D600	WL (nm)		PAT	HLENGTH	
	600		10		
	UNIT		SAM	PLE NAME	
	cells/ml		S-(01	
	FACTOR		MUL	TIPLIER	
	500		100	00	
	OK		Sat	ve	

When select UNIT cells/ml, it will display above parameter setting interface

WL(nm)	Test wavelength
PATHLENGTH	Optical length of cuvette (optical length is less than 10mm)
UNIT	Concentration units
SAMPLE NAME	Sample name
FACTOR	Correction factor
MULTIPLIER	Multiple
OK	Finish parameter setting, back to measure interface automatically
Save	Save parameter (SD card)



S-300 spectrum interface

ABS/%T given wavelength photometry direct-reading measure interface



Wavelegnth scan interface

HOME	I		50	0.0nm	0.000AB	s 12/06/2	8 09:35
ABS/%T	0.01 [:	:	:	:	:	READ
ML SCAN							
	ABS						SET
KINETICS							ZERO
CONC							TOPT
	-0.01 [51		· ·	<u>.</u>	<u>)</u> .	 1000nm	
As As							PNR



S-300 spectrum interface KINETICS dynamics measure interface



CONC concentration regression and measure interface

HOME	500.0nm 0.000ABS 12/06/28	09:35
ABS/%T WL SCAN	$\frac{\overrightarrow{WL}(nm) ABS}{440}$	READ SET
KINETICS	CONCENTRATION	ZERO
CONC		
	SAMP NAME	LOAD
41-41-		PNR



Instructions for installation and use

Appliance / Affair: Spectrophotometer

Accessories installation for S-200/S-220/S-300



CONTENTS

Thermostatic Cell Holder (Part No.: 000A-XB10-0000)	1
Auto Sipper Cell Holder (Part No.: 000B-XB10-0000)	9
Automatic 5-Cell Changer (Part No.: 000C-XB10-0000) 1	19
Micro Cell Holder (100μL) (Part No.: 000D -XB10-0000)2	25
Long Path Cell Holder (100mm) (Part No.: 000F-XB10-0000)2	28
Flow Cell Holder (Part No.: 000G-XB10-0000)	31
Test Tube Holder (Part No.: 000E-VIS20-0000)	36

The instrument can automatically recognize three accessories with the auto sample sipper, auto 5-cell holder, thermostatic cell holder (Only for S-220/S-300).

Remark: When installing accessories, the power of the instrument must turn off.

Thermostatic Cell Holder (Part No.: 000A-XB10-0000)

Overview

The thermostatic cell holder is used to control and maintain the temperature of the sample to be measured. It is usually used for incubate or temperature sensitive sample.



Thermostatic Cell Holder

Description

The thermostatic cell holder is composed of a peltier heating element, a heat sink for cooling of the heating element, a fan for ventilation and a holder for 10mm path length cell. The power is supplied from the main unit DC24V output. The temperature setting is controlled by main unit on-board software or UV Detective PC software (if connected). Water circulation is not required.

Applicable models

The accessory can be used in spectrophotometers model S-220/S-300.

Application

The accessory can be used to measure the transmittance or absorbance of a sample at different temperatures. It can also be used to monitor the change of transmittance or absorbance over time at a designated temperature. Typical application is to measure the enzymatic activity at body temperature 37°C.

Specifications

Thermostatic Cell Holder



Item	Specifications	Remark
Temperature setting range	20- 40°C	
Temperature setting interval	0.1°C	
Temperature accuracy	±0.1°C	
Applicable cell	10mm pathlength cell	Any cell with external dimension of
		12.5 x 12.5mm can be used.

Operating Environment

The thermostatic cell holder generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-220/S-300 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation

• Take out the thermostatic cell from the package.



• Make sure the main power of the instrument is off.





• Dismantle the back cover with screwdriver



• Take the socket of the fan.



• Connect it to the socket of the main unit.





• Insert the fan into the main unit





• Tighten the screws.



• Open the sample compartment cover and remove the screws of the standard cell holder.





• Remove the standard cell holder



• Remove the screw of the bottom plate



• Remove the bottom plate





• Pass the power and data cable through the hole.





 Place the thermostatic cell holder into the sample compartment.





- Locate and align the holes for fixing screws.
- Tighten the fixing screws



 Close the sample compartment. Connect the data cable to the accessory port and the power cable to the DC24V socket at the back of the instrument. Tighten the screw of the socket



 Power on the instrument. The actual temperature is now displaying at the top-left of the home screen.

Temperature





Operation

The thermostatic cell holder can be controlled by the instrument on-board user interface or UV Detective. You can turn on or off of the temperature control and set the temperature.

On-board Interface

Turn on the instrument and you can see the actual cell holder temperature showing on the top-left of the screen. Touch the temperature will enter the temperature setting interface. Default setting of the Temp Control is off. To set the temperature, please touch Temp Control on the right side to turn it on. Then touch the Temperature on the left side and then enter the temperature via the virtual keyboard. Touch Enter and the setting will be effective. You can see the actual temperature is changing to your setting temperature. When the actual temperature reaches to the setting temperature, place the cell containing sample and you can start to measure as normal measurement.

• Touch the Temperature of at the top-left screen and enter the temperature setting screen.

Touch this



 Turn the temperature on and set the temperature (20-40[°]C)



TEMPERATURE (℃): TEMP CONTROL: Setting control temperature Turn on or off temperature control function



Auto Sipper Cell Holder (Part No.: 000B-XB10-0000)

Overview

The auto sipper is used to sip sample solution from a test tube and the sample will flow through the flow cell. The measurement will start automatically. It is usually used for easy/fast sample analysis from any container/sources.



Auto Sipper

Description

The auto sipper is composed of a sipper pump unit and a flow cell. The power is supplied from the main unit DC24V output. The sipper setting is controlled by main unit on-board software or UV Detective PC software (if connected). The sipping actions (Sip, Wash and Return) are performed on the control panel on the sipping pump.

Applicable models

The accessory can be used in spectrophotometers model S-220/S-300 and S-200.

Application

The accessory can be used to measure transmittance or absorbance of a sample from a test tube, a beaker or other container. It can also be used for remote sampling. Typical application is to measure large amount (< 200) of environmental or water samples.

Specifications

Auto Sipper			
Item	Specifications	Remark	
Flow cell volume	150µL		
Flow cell pathlength	10mm		



Sip Time setting range	1-60sec
Delay time setting range	0-200 sec
Purge time setting range	0-60 sec

Operating Environment

The auto sipper generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-220/S-300/S-200 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation



Make sure the main power of the instrument is off.





• Open the sample compartment cover and remove the screws of the standard cell holder.



Remove the standard cell
 holder



 Pass the hose from flow cell through the hole of the bottom plate



Hole





• Place the flow cell into the sample compartment.







• Tighten the fixing screws



- Pull out the black rubber plug of sample compartment and then insert the rubber plug for auto sipper
 - Original rubber



Pull out the original rubber

Insert the rubber plug for auto





• The sipper hose passes the black rubber plug and close the sample compartment cover

Rubber





• Connect the sipper hose to the pump passing the hole.



• Insert the hose to the joint.



• Tighten the nut by a wrench





• Pull the plate towards panel until a rattling sound is heard.



- Push the button to release plate if the auto sipper does not work. Repeat above step again.
- Connect the power cable and data cable to the sipper pump. Tighten the socket to the port.





 Close the sample compartment. Connect the data cable to the accessory port and the power cable to the DC24V socket at the back of the instrument. Tighten the screw of the socket.





 Power on the instrument. The sipper button is now displaying at the top-left of the home screen.
 Sipper



Operation

The auto sipper setting can input from the instrument on-board user interface or UV Detective. You can set the sip time, delay time and the purge time.

On-board Interface

Turn on the instrument and you can see the Sipper button showing on the top-left of the screen. Touch the Sipper button will enter the sipper setting interface. Default setting of the Sip Time is 20 seconds. To set the sip time, touch the Sip Time button and then enter the time via the virtual keyboard. To set delay time, touch the Delay Time button and then enter the delay time via the virtual keyboard. To set the purge time, touch the Purge Time button then enter the purge time via the virtual keyboard. To uch the Home button and continue with other settings.

• Touch the Sipper button at the top-left screen and enter the sipper setting screen.

Touch this




• Touch the Sip Time, Delay Time or Purge Time to set the time for sipper.

HOME	🔮 <mark>Sipper</mark> 500.0nm	0.000ABS 12/06/28 09:35
ABS/%T	SIP TIME(s)	DELAY TIME(s)
	5	2
WL SCAN		
TIME SCAN	0	
CONC		
CONC		
MANAGER		
SETTING		
MANAGER SETTING	ОК	

- SIP TIME (s): sample sipper time (Unit: s)
- DELAY TIME (s):
- delay time between sample sipper and measure
- PURGE TIME (s):
 -): purge air after finish measurement.

This is used to separate the two samples in the tubing. In order to recycle sample or avoid cross contamination.



Procedure

- 1. Touch the "ZERO" button to perform set zero of the instrument.
- 2. Place the sipper hose to the test tube (or other sample container).
- 3. Press "SIPPING" button on the control panel of the sipper pump.
- 4. The sample solution will be sip into the flow cell.
- 5. Press the "READ" button and the measurement will take place
- 6. After measurement, place the sipper hose to the washing solution and the press "WASH" button.
- 7. Repeat above steps for others samples.







Automatic 5-Cell Changer (Part No.: 000C-XB10-0000)

Overview

The Automatic 5-Cell holder is used to hold up to 5 cells with sample and measure each by sequence. It is usually used for fast detection for batch samples.



Automatic 5-Cell Changer

Description

The automatic 5-cell change is composed 5 cell holder arranged in circular manner. A motor is built with the accessory for automatic changing of the cell holder to measure position. Each cell holder can hold one 10mm path length cell. The power is supplied from the main unit DC24V output. The cell holder changing is controlled by main unit on-board software or UV Detective PC software (if connected).

Applicable models

The accessory can be used in spectrophotometers model S-220/S-300 and S-200.

Application

The accessory can be used to measure the transmittance or absorbance from a batch of maximum 5 samples without the need to change the cell by hand. It can also be used to monitor the change of transmittance or absorbance over time from multiple samples such as slow chemical reaction at different initial concentration. Typical application is to place the blank and 3 standard solutions in the cell changer for standardization and quantitation calculation.

Specifications

Automatic 5- Cell Changer



Item	Specifications	Remark
No. of cell holder	5	
Cell changing mechanism	Automatic by software	
Applicable cell	10mm pathlength cell	Any cell with external dimension of
		12.5 x 12.5mm can be used.
		Micro cell is not supported

The automatic 5-cell changer generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-220/S-300/S-200 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation

•



• Make sure the main power of the instrument is off.





• Open the sample compartment cover and remove the screws of the standard cell holder.



 Remove the standard cell holder



• Remove the screw of the bottom plate



• Remove the bottom plate





• Pass the power and data cable through the hole.



 Place the automatic 5-cell holder into the sample compartment.







• Locate and align the holes for fixing screws.

Holes for

fixing



• Tighten the fixing screws



- Close the sample compartment. Connect the data cable to the accessory port and the power cable to the DC24V socket at the back of the instrument. Tighten the screw of the socket
- Power on the instrument. The 5-cell changer button is now displaying at the top-left of the home screen and showing the current position.





5-Cell

Operation

The automatic 5-cell changer can be controlled by the instrument on-board user interface or UV



Detective. You can change the current cell holder position.

On-board Interface

Turn on the instrument and you can see the 5-cell button and the current cell holder position showing on the top-left of the screen. Touch the 5-cell button will enter the cell holder position setting interface. Default setting is position #1. To change the cell position, touch the GOTO CELL button and enter the cell holder position via the virtual keyboard. Place the cell containing sample and you can start to measure as normal measurement.



 Touch GOTO CELL to change the current cell holder position



CELL NUM SETUP: GOTO CELL Cuvette position to be measured (1 to 5) Select specific cuvette position



Micro Cell Holder (100µL) (Part No.: 000D -XB10-0000)

Micro Cell Holder (50μL) (Part No.: 000H -XB10-0000)

Overview

The Micro cell holder (100μ L and 50μ L) is used to hold micro cell of volume down to 100μ L or 50μ L. It is usually used for measuring trace amount of sample or precious sample.



100µL micro cell

Micro Cell Holder

Description

The micro cell holder is designed to use micro cell of sample volume down to 100μ L or 50μ L.

Applicable models

The accessory can be used in spectrophotometers model S-220/S-300 and S-200.

Application

The accessory can be used to measure the transmittance or absorbance from a sample of trace volume of precious sample. Typical application is to measure the concentration of DNA or RNA samples which are very precious and thus only trace amount of sample is available for measurement.

Specifications

Micro Cell Holder		
Item	Specifications	Remark
Micro sample volume	100μL or 50μL	
Beam height	15mm	
Applicable cell	100μL or 50μL micro cell	



The micro cell holder generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-220/S-300/S-200 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation

• Take out the micro cell holder from the package

Micro cell holder



• Open the sample compartment cover and remove the screws of the standard cell holder.



 Remove the standard cell holder





• Place the micro cell holder into the sample compartment.



 Locate and align the holes for fixing screws. Tighten the fixing screws



Holes for fixing

Operation

Place a micro cell containing the sample into the micro cell holder. Close the sample compartment lid and measure the transmittance or absorbance from the on-board software of UV Detective software.



Long Path Cell Holder (100mm) (Part No.: 000F-XB10-0000)

Overview

The long path cell holder is used to hold various long path length cell from 10mm up to 100mm cell. It is usually used for very dilute samples.



Long Path Cell Holder

Description

The micro cell holder is designed to use cell of path length ranges from 10mm, 20mm, 30mm, 40mm, 50mm and 100mm.

Applicable models

The accessory can be used in spectrophotometers model S-220/S-300 and S-200.

Application

The accessory can be used to measure the transmittance or absorbance from a sample of dilute concentration of the analyte. According to Beer's Law, the absorbance of a sample can be increased if the path length is increased. This is useful when the analyte concentration is too low to detect at 10mm path length. Typical application is to measure trace of contaminant such as silica in purified water or environmental samples.

Specifications

Long Path Cell Holder		
Item	Specifications	Remark
Path length	10, 20, 30, 40, 50, 100mm	
Applicable cell	Cell of path length 10, 20,	
	30, 40, 50, 100mm	



The long path cell holder generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-220/S-300/S-200 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation

• Take out the long path cell holder from the package

Long path cell . holder



 Open the sample compartment cover and remove the screws of the standard cell holder.



 Remove the standard cell holder





• Place the long path cell holder into the sample compartment.



 Locate and align the holes for fixing screws. Tighten the fixing screws



Holes for fixing

Operation

Place a long path cell containing the sample into the long path cell holder. Close the sample compartment lid and measure the transmittance or absorbance from the on-board software of UV Detective software.

Adjust the spring position to fit different path length cell
Take out the spring. Insert into the proper position



Flow Cell Holder (Part No.: 000G-XB10-0000)

Overview

The flow cell holder path cell holder is used to allow sample solution to flow from though and measure continuously.



Long Path Cell Holder

Description

The flow cell holder is composed of a flow though cell with volume 150μ L, inlet and outlet tube. External pump such as peristaltic pump is needed.

Applicable models

The accessory can be used in spectrophotometers model S-220/S-300 and S-200.

Application

The accessory can be used to measure the transmittance or absorbance from a sample solution flowing through the flow cell. This is usually use for on-line measurement or remote measurement of a sample not suitable to use with cuvette. Typical application is to measure change of transmittance for fraction collection purpose such as liquid chromatography or continuous flow centrifugation.

Specifications

Long Path Cell Holder		
ltem	Specifications	Remark
Flow cell volume	150µL	
Wavelength range	200-900nm	



Cell material	Quartz
Max. flow rate	10mL/min

The long path cell holder generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-220/S-300/S-200 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation

• Take out the flow cell holder from the package

Flow cell



 Open the sample compartment cover and remove the screws of the standard cell holder.





 Remove the standard cell holder



 Pass the hose from flow cell through the hole of the bottom plate



Hole





• Place the flow cell into the sample compartment.



• Locate and align the holes for the fixing screws



Holes for fixing

• Tighten the fixing screws



 Pull out the black rubber plug Original rubber of sample compartment and then insert the rubber plug for auto sipper

r



For auto sipper

Pull out the original rubber Insert the rubber plug for auto





• The sipper hose passes the black rubber plug and close the sample compartment cover

Rubber



Operation

Connect the inlet tube to a peristaltic pump or other external pump. Feed the sample solution to the flow cell. Measure the transmittance or absorbance from the on-board software of UV Detective software.



Test Tube Holder (Part No.: 000E-VIS20-0000)

Overview

The test tube holder is used for the direct measurement of samples in a test tube with no need to transfer to a cuvette.



Test tube holder

Description

The test tube holder is composed of a spring mechanism which can adjusts and accepts test tubes with diameters from 9-22mm automatically. And it has a high ceiling cover accommodates even the tallest test tube.

Applicable models

The accessory can be used in spectrophotometers model S-200.

Application

The accessory can be used to measure the transmittance or absorbance directly from a test tube, which is very convenient.

Specifications			
Test Tube Holder			
Item	Specifications	Remark	
Tube size			
Diameter	9-22mm		
Length	70-150		
Beam height	21.5mm		

Specifications



The test tube holder generally operates under the same environmental conditions as those specified in the manual of the Boeckel + Co (GmbH + Co) KG S-200 Spectrophotometer.

Unpacking

Unpack the box and take out the items inside. Check all items are available.

Installation

 Take out the test tube holder and the accessory from the package



• Open the sample compartment cover and remove the screws of the standard cell holder.







 Put the test tube holder in the sample chamber and flat it by inserting the two position pins in the holes.



Load and tighten the fixing screws.



• Remove the turn door using the allen wrench in the accessory.







 Install the turn door in the accessory. Load the fixing screw and push back the turn door, then you can tighten the fixing screw by allen wrench.







Operation

Place a test tube containing the sample into the test tube holder directly. Close the sample compartment lid and measure the transmittance or absorbance from the on-board software of UV Detective software.