

BOECO DESKTOP INCUBATOR I-18 AND UNIVERSAL DESKTOP OVEN U-18

The I-18/U-18 are featured with unique design, modern appearance, and advanced manufacturing techniques. They have an accurately temperature control system, daily maintains and operates easily. Those Incubator/Oven are widely used in bacteria & microbial scientific research, medical science, pharmacy, biological and biochemical industries etc

Features:

- ▶ Vertical design, box body made by high quality cold-rolled steel, the surface is treated by electrostatic spraying process, stable and beautiful appearance.
- ▶ LCD fuzzy-PID control temperature controller with overtemperature protection
- ▶ Stainless Steel inner chamber, Circular arc transition corners, adjustable shelf, and easy cleaning.
- ▶ The external chamber is processed by cold rolling steel plate after anti-corrosive coating.
- ▶ Door with double toughened glass window enables large angle observation. Elastic adjustable lock type handle. High temperature silica gel sealing strip, to ensure good sealing and prevent heating loss.



Specification I-18:

Volume:	18 l
Heating power:	200W
Operating Temp. Range:	Room Temperature + 5°C... + 65°C
Temp. Uniformity:	± 0,5
Temp. Fluctuation:	± 0,5
Temperature Setting:	0,1°C
Time range:	1 -9999min, continuous.
Max. Number of shelves:	2
Shelf Load:	15 kg
Internal Dimension:	260 x 260 x 260 mm (L x W x H)
External Dimension:	460 x 510 x 695 mm (L x W x H)
Weight net:	20 kg

Code	Description
BOE 8042000	BOECO Desktop Incubator I-18, 220V, 50-60 Hz incl. 2 stainless steel shelves
BOE 8042100	BOECO Desktop Incubator I-18, 115V, 60 Hz incl. 2 stainless steel shelves

Specification U-18:

Volume:	18 l
Heating power:	600W
Operating Temp. Range:	Room Temperature + 10°C... + 300°C
Temp. Uniformity:	± 0,5
Temp. Fluctuation:	± 0,5
Temperature Setting:	0,1°C
Time range:	1 -9999 min, continuous.
Max. Number of shelves:	2
Shelf Load:	15 kg
Internal Dimension:	260 x 260 x 260 mm (L x W x H)
External Dimension:	460 x 510 x 695 mm (L x W x H)
Weight net:	20 kg

Code	Description
BOE 8043000	BOECO Universal Desktop oven U-18, 220V, 50-60 Hz incl. 2 stainless steel shelves
BOE 8043100	BOECO Universal Desktop oven U-18, 115V, 60 Hz incl. 2 stainless steel shelves