

PIPETTE AIDS



**BOECO PIPETTE
FILLER BALLS**

BOECO PIPETTE FILLER BALLS

natural rubber, red, with 3 valves

Cat. No.	Description
BOE 100	Model STANDARD, short cone for pipettes up to 20 ml
BOE 101	Model UNIVERSAL, long cone for all pipettes



**PIPETTE PUMPS,
PP**

PIPETTE PUMPS, PP

can be operated with one hand. Glass and plastic pipettes fit easily. Precise filling or emptying is done by rotating their actuator-wheels forward or backward. For rapid emptying simply depress plunger.

Cat. No.	Description
BOE 120	0 - 2 ml, blue
BOE 121	0 - 10 ml, green
BOE 122	0 - 25 ml, red



**BOECO PIPETTE
CONTROLLER**

BOECO PIPETTE CONTROLLER

offers safer, more convenient filling and dispensing, with glass and plastic pipettes. With one sensitive lever filling and delivery are easily controlled. Pipetting with ease and comfort.

- » Easy and tireless handling for left and right-handed people
- » An especially designed intake cone allows the use of pipettes from 0,1 ml to 100 ml
- » The perfectly harmonized valve unit ensures effortless liquid intake
- » Precise and high sensitive set of the meniscus
- » An easy replaceable hydrophobic membrane filter protects the interior from penetrating liquids
- » Completely autoclavable at 121°C (2 bar)

Cat. No.	Description
BOE 160	BOECO Manual Pipette controller, complete
BOE 1670650	Membrane filter, 3 µm, sterile, hydrophobic, PTFE, pack of 10 pcs.
BOE 1670211	Suction bellows / Intake element
BOE 1670202	Intake cone, silicone
BOE 1655590	Valve unit