

STEREO MICROSCOPES

BOECO ZOOM STEREO MICROSCOPE MODEL BSZ-405

Specification BSZ-405:

Optical system:	Binocular head inclined at 45° Interpupillary Distance 55 - 75 mm Eyepiece dioptic adjustment range $\pm 5D$
Working distance:	97 mm
Field of View:	20,0 to 4.4 mm diam.
Eyepiece:	EWF 10x/20 mm eyepieces
Zoom Optics:	1X to 4,5x , Zoom ratio 1:4,5
Illumination:	2-way illumination LED 1,2W Reflected (upper) illumination LED 0,7W Transmitted (lower) illumination
Stage Plate:	Frosted glass stage plate 90 mm and black/white reversible stage plate
Base:	Stable 262 x 185 mm base, fitted with anti-skid rubber feet. Main Body 360° rotatable
Power:	100-240 V, 50/60 Hz
Dimensions:	185 (W) x 262 (L) x 315 (H) mm, Net weight 4,3 / 4,4 kg

Code	Description	Head	Plug
BOE 4050.000	BSZ-405	Binocular	Euro
BOE 4050.001	BSZ-405	Binocular	US
BOE 4050.300	BSZ-405/T	Trinocular	Euro
BOE 4050.301	BSZ-405/T	Trinocular	US

optional

Code	Description
BOE 3460.020	Wide Field Eyepiece WF 20x/12 with diopter adjustment
BOE 4050.002	Auxiliary Objective 2x, for BSZ-405 and BS-80
BOE 4050.007	Auxiliary Objective 0,75x, for BSZ-405 and BS-80

BOECO ZOOM STEREO MICROSCOPE MODEL BST-606

Specification BST-606:

Optical system:	Binocular head inclined at 45° Interpupillary distance 55 - 75 mm Eyepiece dioptic adjustment range $\pm 5D$
Working distance:	115 mm
Field of View:	28,0 to 4,5 mm diam.
Eyepiece:	EWF 10x/22 mm eyepieces with eyeshadow
Zoom Optics:	0,8x to 5x, Zoom ratio 6,3:1
Illumination:	2-way illumination switch allows upper and lower illumination Transmitted (lower) LED light source Incident (upper) LED light source
Stage Plate:	100 mm $\varnothing$ frosted glass stage plate and black/white reversible stage plate
Base:	Stable 300x225 mm base, fitted with anti-skid rubber feet. Main Body 360° rotatable
Power:	100-240 V, 50/60 Hz
Dimensions:	300 (W) x 225 (L) x 300 (H) mm, Net weight 5,4 kg

Code	Description	Head	Plug
BOE 6060.000	BST-606	Binocular	EURO
BOE 6060.001	BST-606	Binocular	US
BOE 6060.300	BST-606/T	Trinocular	Euro
BOE 6060.301	BST-606/T	Trinocular	US

optional

BOE 3460.020	Wide Field Eyepiece WF 20x/12 w/ diopter adjustment
BOE 6060.002	Auxiliary Objective 2x, for BST-606
BOE 6060.005	Auxiliary Objective 0,5x, for BST-606



BSZ-405 TRINOCULAR WITH DIGITAL MICROSCOPE CAMERA B-CAM10

